

	REVIEWER	VIEWER	TECHNICIAN	USER	OWNER
Note	Available only on P/E/T		Available only on P/E/T		
Workspace	N/A	No abilities	N/A	Create & edit tags available in workspace Create projects	<b>Everything in User plus:</b> Create & edit workspaces and all related settings (tags, automations, groups) View audit trails
Projects, Experiments, Tasks (in short "P/E/T")	View P/E/T Comment task Comment task protocol steps & results Create task signature requests & sign task	View P/E/T	<b>Everything in Viewer plus:</b> Comment projects Comment task Change task status Edit task assigned items & consume item stock Complete task protocol steps & checklist items Comment task protocol steps & results Edit form fields on task protocol step & submit form Create task signature requests & sign task	<b>Everything in Technician plus:</b> Create & edit P/E/T Edit task protocol Create & edit task results Add form to task protocol step & Unlock submitted form	<b>Everything in User plus:</b> Edit access to P/E/T Edit task assigned user
Project insights	N/A	View insights for Projects with access	N/A	View insights for Projects with access	View insights for Projects with access
Inventories	N/A	View	N/A	<b>Everything in Viewer plus:</b> Create & edit items Archive items Create & edit inventory filters	<b>Everything in User plus:</b> Edit access Share Create & edit inventory columns
Item locations	N/A	View	N/A	<b>Everything in Viewer plus:</b> Create & edit	<i>Same as User</i>
Forms	N/A	View	N/A	<b>Everything in Viewer plus:</b> Create & edit form draft	<b>Everything in User plus:</b> Edit access Publish Archive
Protocol templates	N/A	View	N/A	<b>Everything in Viewer plus:</b> Create & edit protocol template draft	<b>Everything in User plus:</b> Edit access Publish Archive

<b>Label templates</b>	<i>N/A</i>	View	<i>N/A</i>	<b>Everything in Viewer plus:</b> Create & edit	<i>Same as User</i>
<b>Reports</b>	<i>N/A</i>	View	<i>N/A</i>	<b>Everything in Viewer plus:</b> Create & edit	<i>Same as User</i>

	<b>Org Admin abilities</b>
<b>Organization administrator</b>	<ul style="list-style-type: none"> <li>Manage organization members</li> <li>Manage organization settings</li> <li>Access system log</li> <li>View &amp; create &amp; invite oneself to any workspace</li> <li>Manage integrations</li> <li>Revoke &amp; delete any task signature request</li> </ul>