

	Reviewer Only on P/E/T	Viewer	Technician Only on P/E/T	User	Owner
Workspace	N/A	No abilities	N/A	<ul style="list-style-type: none"> - Create & edit tags available in workspace - Create projects 	Everything in User plus: <ul style="list-style-type: none"> - Create & edit workspaces and all related settings (tags, automations, groups) - View audit trails
Projects, Experiments, Tasks (in short "P/E/T")	<ul style="list-style-type: none"> - View P/E/T - Comment task - Comment task protocol steps & results - Create task signature requests & sign task 	View P/E/T	Everything in Viewer plus: <ul style="list-style-type: none"> - Comment projects - Comment task - Change task status - Edit task assigned items & consume item stock - Complete task protocol steps & checklist items - Comment task protocol steps & results - Edit form fields on task protocol step & submit form - Create task signature requests & sign task 	Everything in Technician plus: <ul style="list-style-type: none"> - Create & edit P/E/T - Edit task protocol - Create & edit task results - Add form to task protocol step & Unlock submitted form 	Everything in User plus: <ul style="list-style-type: none"> - Edit access to P/E/T - Edit task assigned user
Project insights	N/A	View insights for Projects with access	N/A	View insights for Projects with access	View insights for Projects with access
Inventories	N/A	View	N/A	Everything in Viewer plus: <ul style="list-style-type: none"> - Create & edit items - Archive items - Create & edit inventory filters 	Everything in User plus: <ul style="list-style-type: none"> - Edit access - Share - Create & edit inventory columns
Item locations	N/A	View	N/A	Everything in Viewer plus: <ul style="list-style-type: none"> - Create & edit locations 	Same as user
Forms	N/A	View	N/A	Everything in Viewer plus: <ul style="list-style-type: none"> - Create & edit form drafts 	Everything in User plus: <ul style="list-style-type: none"> - Edit access - Publish forms - Archive forms
Protocol templates	N/A	View	N/A	Everything in Viewer plus: <ul style="list-style-type: none"> - Create & edit protocol template draft 	Everything in User plus: <ul style="list-style-type: none"> - Edit access - Publish templates - Archive templates
Label templates	N/A	View	N/A	Everything in Viewer plus: <ul style="list-style-type: none"> - Create & edit labels 	Same as user
Reports	N/A	View	N/A	Everything in Viewer plus: <ul style="list-style-type: none"> - Create & edit labels 	Same as user

Organization admin abilities

- Manage organization members
- Manage organization settings
- Access system log
- View & create & invite oneself to any workspace
- Manage integrations
- Revoke & delete any task signature request